

MANUAL

Scythe can be run minimally with:

```
scythe -a adapter_file.fasta -o trimmed_sequences.fasta sequences.fastq
```

Run scythe the adaptor file the outfile the infile

By default, Illumina's quality scheme (pipeline > 1.3) is used. Sanger or Solexa (pipeline < 1.3) qualities can be specified with `-q`:

```
scythe -a adapter_file.fasta -q solexa -o trimmed_sequences.fasta sequences.fastq
```

Q2: Based on the the above – indicate which words go in which box.

/scythe -a -q -o

- a) `$dir_scy` (variable with the path to the software scythe)
- b) `1_scy_human001_pe1.fastq` (outfile name)
- c) `adaptors.fa`
- d) `human001_pe1.fastq`
- e) `sanger`

MANUAL

Scythe can be run minimally with:

```
scythe -a adapter_file.fasta -o trimmed_sequences.fasta sequences.fastq
```

Run scythe the adaptor file the outfile the infile

By default, Illumina's quality scheme (pipeline > 1.3) is used. Sanger or Solexa (pipeline < 1.3) qualities can be specified with `-q`:

```
scythe -a adapter_file.fasta -q solexa -o trimmed_sequences.fasta sequences.fastq
```

Q2: Based on the the above – indicate which words go in which box.

```
$dir_scy/scythe -a adaptors.fa -q sanger -o 1_scy_human001_pe1.fastq human001_pe1.fastq
```

- a) `$dir_scy` (variable with the path to the software scythe)
- b) `1_scy_human001_pe1.fastq` (outfile name)
- c) `adaptors.fa`
- d) `human001_pe1.fastq`
- e) `sanger`

MANUAL

Scythe can be run minimally with:

```
scythe -a adapter_file.fasta -o trimmed_sequences.fasta sequences.fastq
```

Run scythe the adaptor file the outfile the infile

By default, Illumina's quality scheme (pipeline > 1.3) is used. Sanger or Solexa (pipeline < 1.3) qualities can be specified with `-q`:

```
scythe -a adapter_file.fasta -q solexa -o trimmed_sequences.fasta sequences.fastq
```

Q2: Based on the the above – indicate which words go in which box.

/scythe -a -q -o

- a) `$dir_scy` (variable with the path to the software scythe)
- b) `1_scy_human001_pe2.fastq` (outfile name)
- c) `adaptors.fa`
- d) `human001_pe2.fastq`
- e) `sanger`

MANUAL

Scythe can be run minimally with:

```
scythe -a adapter_file.fasta -o trimmed_sequences.fasta sequences.fastq
```

Run scythe the adaptor file the outfile the infile

By default, Illumina's quality scheme (pipeline > 1.3) is used. Sanger or Solexa (pipeline < 1.3) qualities can be specified with `-q`:

```
scythe -a adapter_file.fasta -q solexa -o trimmed_sequences.fasta sequences.fastq
```


Q2: Based on the the above – indicate which words go in which box.

```
$dir_scy/scythe -a adaptors.fa -q sanger -o 1_scy_human001_pe2.fastq human001_pe2.fastq
```

- a) `$dir_scy` (variable with the path to the software scythe)
- b) `1_scy_human001_pe2.fastq` (outfile name)
- c) `adaptors.fa`
- d) `human001_pe2.fastq`
- e) `sanger`

MANUAL

Run sickle, Paired end mode, input forward (pe1) reads, input reverse (pe2) reads



Q3: Based on the the above – indicate which words go in which box. We want to trim to a quality of 20 and length of 20

sickle -f -r -t -o -p -s -q -l

- a) sanger
- b) 1_scy_human001_pe1.fastq
- c) 1_scy_human001_pe2.fastq
- d) \$dir_sic
- e) 2_sic_human001_pe1.fastq
- f) 2_sic_human001_pe2.fastq
- g) 2_sic_human001_sing.fastq
- h) 20
- i) pe
- j) 20

MANUAL

Run sickle, Paired end mode, input forward (pe1) reads, input reverse (pe2) reads

```
sickle pe -f input_file1.fastq -r input_file2.fastq -t sanger \
-o trimmed_output_file1.fastq -p trimmed_output_file2.fastq \
-s trimmed_singles_file.fastq -q 12 -l 15
```

Annotations:

- pe: Paired end mode
- f input_file1.fastq: input forward (pe1) reads
- r input_file2.fastq: input reverse (pe2) reads
- t sanger: Data files = sanger format
- o trimmed_output_file1.fastq: forward pe1outfile
- p trimmed_output_file2.fastq: reverse pe2outfile
- s trimmed_singles_file.fastq: outfile_singletons
- q 12 -l 15: Trim if quality is <12, remove if read < 15 bases

Q3: Based on the the above – indicate which words go in which box. We want to trim to a quality of 20 and length of 20

```
$dir_sic/sickle pe -f 1_scy_human001_pe1.fastq -r 1_scy_human001_pe1.fastq -t sanger \
-o 2_sic_human001_pe1.fastq -p 2_sic_human001_pe2.fastq \
-s 2_sic_human001_sing.fastq -q 20 -l 20
```

Two steps – 1 (has to be done for pe1 & pe2)

MANUAL

Run BWA, align mode, genome reference, input file output file

The diagram shows the command `bwa aln ref.fa short_read.fq > aln_sa.sai`. Red arrows point from the labels above to the components: 'bwa' from 'Run BWA', 'aln' from 'align mode', 'ref.fa' from 'genome reference', 'short_read.fq' from 'input file', and 'aln_sa.sai' from 'output file'.

Based on the the above indicate which words go in which box

/bwa dir_ref/ >

aln
 2_sic_human001_pe1.fastq
 dir_bwa
 3_bwa_human001_pe1.sai
 chr14.fa

Two steps – 1 (has to be done for pe1 & pe2)

MANUAL

Run BWA, align mode, genome reference, input file output file

The diagram shows the command line `bwa aln ref.fa short_read.fq > aln_sa.sai`. Red arrows point from the labels above to the components: 'bwa' to 'bwa', 'aln' to 'aln', 'ref.fa' to 'ref.fa', 'short_read.fq' to 'short_read.fq', and 'aln_sa.sai' to 'aln_sa.sai'.

Based on the the above indicate which words go in which box

dir_bwa/bwa aln dir_ref/chr14.fa 2_sic_human001_pe1.fastq >
3_bwa_human001_pe1.sai

aln
2_sic_human001_pe1.fastq
dir_bwa
3_bwa_human001_pe1.sai
chr14.fa

Two steps – 1 (has to be done for pe2 now)

MANUAL

Run BWA, align mode, genome reference, input file output file

The diagram shows the command line: `bwa aln ref.fa short_read.fq > aln_sa.sai`. Red arrows point from the labels above to the components: 'bwa' (Run BWA), 'aln' (align mode), 'ref.fa' (genome reference), 'short_read.fq' (input file), and 'aln_sa.sai' (output file).

Based on the the above indicate which words go in which box

/bwa
 dir_ref/

 >

aln
 2_sic_human001_pe2.fastq
 \$dir_bwa
 3_bwa_human001_pe2.sai
 chr14.fa

Two steps – 1 (has to be done for pe1 & pe2)

MANUAL

Run BWA, align mode, genome reference, input file output file

The diagram shows the command `bwa aln ref.fa short_read.fq > aln_sa.sai`. Red arrows point from the labels above to the components: 'bwa' (Run BWA), 'aln' (align mode), 'ref.fa' (genome reference), 'short_read.fq' (input file), and 'aln_sa.sai' (output file).

Based on the the above indicate which words go in which box

`$dir_bwa/bwa` `aln` `dir_ref/chr14.fa` `2_sic_human001_pe2.fastq` `>`
`3_bwa_human001_pe2.sai`

aln
 2_sic_human001_pe1.fastq
 \$dir_bwa
 3_bwa_human001_pe1.sai
 chr14.fa

Two steps -2

MANUAL

run BWA, paired end mode, reference,

INPUT FILES

- sai files from step 3
- fastq files from step 2

output

```
bwa sampe ref.fa aln_sa1.sai aln_sa2.sai read1.fq read2.fq > aln-pe.sam
```

Based on the the above indicate which words go in which box

```
$dir bwa/bwa sampe $dir ref/chr14.fa
```

>

2_sic_human001_pe1.fastq

```
2 sic human001 pe2.fastq
```

3 bwa human001 pe1.sai

3 bwa human001 pe2.sai

```
4 bwa human001 pe12.sam
```

Two steps -2

MANUAL

run BWA, paired end mode, reference,

INPUT FILES

- sai files from step 3
- fastq files from step 2

output

```
bwa sampe ref.fa aln_sa1.sai aln_sa2.sai read1.fq read2.fq > aln-pe.sam
```

Based on the the above indicate which words go in which box

```
$dir_bwa/bwa sampe $dir_ref/chr14.fa 3_bwa_human001_pe1.sai
```

```
3_bwa_human001_pe2.sai 2_sic_human001_pe1.fastq
```

```
2_sic_human001_pe2.fastq > 4_bwa_human001_pe12.sam
```

```
2_sic_human001_pe1.fastq
```

```
2_sic_human001_pe2.fastq
```

```
3_bwa_human001_pe1.sai
```

```
3_bwa_human001_pe2.sai
```

```
4_bwa_human001_pe12.sam
```

Sam to Bam file

Lets convert our sam file to a bam file using picard tools

MANUAL

```
java -jar SortSam.jar \  
  I=input.sam \  
  O=output.bam \  
  SO=coordinate
```

We also have to add: CREATE_INDEX=true

```
java -jar /SortSam.jar I=  
O= SO=coordinate CREATE_INDEX=true
```

\$dir_pic

5_bwasort_human001_pe12.bam

4_bwa_human001_pe12.sam

Sam to Bam file

Lets convert our sam file to a bam file using picard tools

MANUAL

```
java -jar SortSam.jar \  
  I=input.sam \  
  O=output.bam \  
  SO=coordinate
```

We also have to add: CREATE_INDEX=true

```
java -jar $dir_pic/SortSam.jar I=4_bwa_pe12.sam  
O=5_bwasort_pe12.bam SO=coordinate CREATE_INDEX=true
```

\$dir_pic

5_bwasort_human001_pe12.bam

4_bwa_human001_pe12.sam

Filtering out poorly mapping reads

CGATACGATATGGTCATGACACGAGTACGATCGATGGTCATGACACGAGTATCGGATGGTCAT
ACGATATGGTCATGACACGAGTACGATCGATGGTCATGACACGAGTATCGGATGGTCAT
TATGGTCATGACACGAGTACGATCGATGGTCATGACACGAGTATCGGATGGTCAT
GTATGACACGAGTACGATCGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
CATGACACGAGTACGATCGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
TGGTAGGATCGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
TGGTATCGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
GGTATCGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
TCGATCGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
ATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
ATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
CCAGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
AGGTGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT
GGTGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGTATCGGATGGTCATGACACGAGT

What if wanted to excluded any unmapped reads or reads with low mapping scores?

Filter out unmapped reads

Filter out reads if mapping quality < 30

`samtools view -F 4 -q 30 -b mybamfile > myoutfile`

`/samtools view -F 4 -q 30`

`-b`

`>`

`5_bwasort_human001_pe12.bam`

`$dir_sam`

`6_bwafil_human001_pe12.bam`

What if wanted to excluded any unmapped reads or reads with low mapping scores?

Filter out unmapped reads

Filter out reads if mapping quality < 30

```
samtools view -F 4 -q 30 -b mybamfile > myoutfile
```

```
$dir_sam/samtools view -F 4 -q 30 -b 5_bwasort_human001_pe12.bam
```

```
> 6_bwafil_human001_pe12
```